

SENSITIVITY SPECIFICITY IN LOCK-IN AMPLIFIERS

by J. L. Scott

Rev. 2 12/93

DL INSTRUMENTS, LLC

233 Cecil A. Malone Dr., Ithaca, NY 14850 Phone 607-277-8498; Fax 607-277-8499

Manufacturers of lock-in amplifiers (LIA) often provide very high sensitivity ranges (e.g., 10 nanovolts full scale) on their instruments and then make this a key sales feature. This is a misleading tactic, since the input-referred spectral noise (e_n) of the LIA will make this range practically useless.

Consider one such specification for a competitive instrument, which lists $e_n = 7 \text{ nV}/\sqrt{\text{Hz}}$ (typical) at 1 kHz and a maximum time constant (T) of 100 seconds with only a 6 dB/octave rolloff available. For this LIA the narrowest possible equivalent noise bandwidth (B), which will yield the least output fluctuation, will be $1/(4T)$ or 0.0025 Hz. The signal fluctuation will be the inverse of the output-referred signal-to-noise ratio. Therefore, per ITHACO Application Note IAN49 "*Speed/Accuracy Tradeoff When Using A Lock-In Amplifier To Measure Signal In The Presence Of Random Noise*", for a full scale input (e_s) of 10 nV on the 10 nV sensitivity range we have a reproducibility given by:

$$\frac{1}{\text{SNR}} = \frac{e_n \sqrt{B}}{e_s} = \frac{7 \text{ nV}}{10 \text{ nV}} \sqrt{0.0025} = 0.035, \text{ rms}$$

This figure (1/SNR) represents the rms output variation expressed as a fraction of the signal reading. It is mathematically identical to the standard deviation (σ) of the measurement. The output peak-to-peak fluctuation will be roughly six times the rms value (6σ). Therefore, for the example above, we would

observe a 21% (2.1 nV) peak-to-peak variation in the measurement ($.035 \times 6 \times 100\%$) with a fluctuation rate on the order of several minutes. On an analog meter, we could read this sort of result just as well on the 100 nV scale (will observe 10% of full scale). For a digital readout the same accuracy can be obtained on the 1 μV sensitivity range due to the higher resolution available.

The point is this: the 10 nV range buys you *nothing*. This is the reason ITHACO does not provide a 10 nV range on its instruments. For example, using a DYNATRAC® instrument with $5 \text{ nV}/\sqrt{\text{Hz}}$ max noise, $T_{\text{max}} = 125$ seconds, 12 dB/octave rolloff ($B = 1/(8T) = .001$ Hz) and 100 nV sensitivity will yield a reproducibility relative to 100 nV of 0.0016, or a 6σ peak-to-peak output fluctuation of 1% of full scale (1.0 nV).

The key feature which places an ultimate limit on resolution is the self noise of the lock-in rather than its output scaling. Even this limit may not be attained in less carefully designed instruments. The noise is usually expressed in volts with the input shorted (e_n). Be cautioned that if you are working with high impedance sources, the effect of LIA input self *current* noise (not always specified) may far exceed that of the voltage noise. You can gauge current noise when comparing lock-ins by connecting a scope directly to the open-circuited signal input.